

Constant Temperature and Humidity Incubator


Features:

- * UV lamp for sterilization.
- * The starting up time can be preset.
- * Universal castors for easy movement.
- * Double door design and magnetic door sealing.
- * Microprocessor temperature control.

Safety Function:

- * Over temperature protection, temperature sensor malfunction alarm.
- * High and low water level alarm, automatic heating off for low water level.

- * LED display
- * LCD display
- * Electrical power outlet in the chamber.
- * CFC-free refrigerant with high efficiency.
- * The working chamber adopts round angle structure, easy to clean.
- * Parameters can be saved automatically when power off or system error.



Humidity area illustration




Micro Printer

Notice: Please don't use in low temperature ($\leq 20^{\circ}\text{C}$) for long time ($\geq 160\text{h}$).

Technical Parameters:

| Model | KLI-160 | KLI-200 | KLI-250 | KLI-300 | KLI-400 |
|----------------------------|--|--------------|--------------|--------------|--------------|
| Capacity | 160L | 200L | 250L | 300L | 400L |
| Temp. Range | 0~65°C | | | | |
| Temp. Fluctuation | $\pm 1^{\circ}\text{C}$ | | | | |
| Temp. Stability/Uniformity | $\pm 2^{\circ}\text{C}$ | | | | |
| Humidity Range | 50~90% RH ($\pm 5\sim 8\%$ RH) | | | | |
| Sterilization Type | UV lamp | | | | |
| Shelves | 3 pcs/adjustable | | | | |
| Consumption | 300W | 350W | 430W | 480W | 750W |
| Power Supply | AC110/220V $\pm 10\%$, 50/60Hz | | | | |
| Standard Accessory | Micro printer(for HTII series), humidifier | | | | |
| Internal Size(W*D*H) mm | 450*420*850 | 450*450*1000 | 480*490*1070 | 520*500*1170 | 580*540*1270 |
| External Size(W*D*H) mm | 590*620*1400 | 590*640*1550 | 620*680*1620 | 660*690*1720 | 720*740*1820 |
| Package Size(W*D*H) mm | 750*780*1600 | 750*800*1740 | 780*840*1810 | 820*850*1910 | 880*900*2010 |
| Gross Weight(kg) | 120 | 138 | 147 | 157 | 170 |



Humidity area illustration