

# DS-FS20E

[desmonscientific.com/](http://desmonscientific.com/)





### BLOOD BANK CONTROLLER FEATURES

End user interface	Display LED
Resolution (°C)	0,1
Maximum Number of Probes	3
Internal Led light	yes
Acoustic alarm	yes
Hot gas advanced defrost	yes

### ALARMS AND FAILURES (ACOUSTIC AND VISUAL)

HACCP function	yes
"Power failure" (back up battery) Alarm	yes
"High / low temperature" Alarm	yes
"Open door" Alarm	yes
"Sensor fault" Alarm	yes
Clogged condenser Alarm	yes
Temperature datalogger with USB port	yes

### CONNECTIONS

RS 485 interface	optional
Porta Ethernet (LAN)	optional

---

## SAFETY

---

Switch on/off protected by password	yes
-------------------------------------	-----

---

Dry Contact for alarms remoting	yes
---------------------------------	-----

---

Freezing and overheating safety device	yes
----------------------------------------	-----

---

## OPTIONALS

---

Rechargeable back-up batteries	optional
--------------------------------	----------

---

## STRUCTURE

---

Exterior dimensions WxDxH (mm)	600x600x1560
--------------------------------	--------------

---

Capacity (L)	200
--------------	-----

---

N of blood bags (450ml)	128
-------------------------	-----

---

Insulation (mm)	60
-----------------	----

---

Temperature range (°C)	+4°C
------------------------	------

---

Max. ambient temperature (°C)	43
-------------------------------	----

---

Type of cooling unit	Ventilated
----------------------	------------

---

Power supply (V)	220V/240V
------------------	-----------

---

Frequency (Hz)	50
----------------	----

---

Power consumption (W)	245
-----------------------	-----

---

Weight (KG)	94
-------------	----

---

Shipping weight (KG)	105
----------------------	-----

---

Shipping dimensions WxDxH (mm)	680x680x1790
--------------------------------	--------------

---

No. of compressors.	1
---------------------	---

---

Doors	1 glass door
-------	--------------

---

Automatic (smart) defrost and evaporation	yes
-------------------------------------------	-----

---

Castors	yes
---------	-----

---

Chart recorder	yes
----------------	-----

---

Drawers (pcs)	4 s/steel drawers
Lock	yes
Interior material	s/steel AISI304
Exterior material	White powder coated steel
Gas	R404a

Price: € 0,00